

NISUS[®] Negative Pressure Wound Therapy System

A powerful yet portable device, allowing patients to achieve results while still maintaining mobility.

- ✓ FDA approved for use in hospitals, long-term care facilities and home-care settings
- ✓ Lightweight and portable, with up to 24 hours of battery life (can extend to 48 hours with sleep mode)
- ✓ User-friendly interactive interface with easy-to-read LCD screen



Nisus NPWT Wound Dressing Kit

Nisus NPWT Pump

Nisus NPWT Canisters

The Cork Medical Nisus Negative Pressure Wound Therapy (NPWT) System is designed to provide gentle powered suction to treat and promote wound healing. The Nisus pump is portable and delivers precise pressure for optimal power and performance to provide a universal solution for all wounds across the patient's full continuum of care.

NISUS NPWT SYSTEM FEATURES

The Cork Medical Nisus Negative Pressure Wound Therapy System is indicated for use in patients who would benefit from negative pressure wound therapy, particularly as the device may promote wound healing by the removal of excess exudates, infectious material and tissue debris.

- **Recommended for patients with:**
 - Chronic, acute, traumatic, subacute and dehisced wounds
 - Partial-thickness burns
 - Pressure ulcers, diabetic ulcers and venous ulcers
 - Flaps and grafts
- **Step-by-step troubleshooting guide for alarms, such as:** Canister full, vacuum leak and blockage
- **IV pole accessible**
- **Powerful suction capacity: ~8-10 L/min**
- **Multiple unique features for clinical enhancement of negative pressure wound therapy, such as:**
 - Patient runtime
 - Sleep feature
 - Alarm log
 - Lock restriction
 - Variable Intermittent and Continuous Therapy Modes
- **Convenient rear compartment battery access and replacement**

COMPONENT SPECS

Composed of three essential components that actively work together to promote wound healing, the Nisus NPWT System includes the Nisus NPWT Pump, the Nisus NPWT Canisters and the Nisus NPWT Wound Dressing Kit.

Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use. Use only Cork Medical components with the Nisus Negative Pressure Wound Therapy System.

- **Battery charger and carrying bag included**
- **Medicare HCPCS Code: E2402**



| Model No. | Component |
|------------------|----------------------------------|
| CMPP-100 | Nisus NPWT Pump |
| CPC-250 | Nisus NPWT 250mL Canister |
| CPC-500 | Nisus NPWT 500mL Canister |
| WDK-250 | NPWT Wound Dressing Kit (Small) |
| WDK-300 | NPWT Wound Dressing Kit (Medium) |
| WDK-400 | NPWT Wound Dressing Kit (Large) |
| CLEANSE-KIT-A-18 | NPWT Cleanse Therapy Kit |
| NIS-BRDR-A-16 | Bridge Foam |
| FOAM-WHITE-A-16 | Optifoam White Foam |
| CSP-100-ES-EXT | 24" Extension Tubing |
| CWY-100 | Y-Connector |

In Partnership with

For more information, contact the InfuSystem NPWT team at 1-844-534-1997 or npwt@infusystem.com.

